



19301

Dextran ex. Leuconostoc Sp. - Tech 500 (450-550)

Part B

CAS : 9004-54-0

Molecular Formula : (C₆H₁₀O₅)_n

Molecular Weight : 5,00,000

Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Powder
Loss on drying	max. 7%
Specific rotation [?]	+195° to +203° (C=2, H ₂ O)
Average Molecular weight	450 - 550 KDa
Not for medicinal use	

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	39139090 (GST 18%)
100 Gms	39139090 (GST 18%)
25 Gms	39139090 (GST 18%)
500 Gms	39139090 (GST 18%)
Type of Packing	
10 Gms	Glass Bottle

Available Packages

10 Gms

25 Gms

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com