



19428

Creatinine extrapure AR, 99%

Part B

CAS : 60-27-5

Molecular Formula : C₄H₇N₃O

Molecular Weight : 113.12

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline Powder
Solubility (Turbidity) 5% aq. solution	Clear
Solubility (Colour) 5% aq. solution	Colourless
Assay (NT) (On dry basis)	min. 99%
pH (5% aq. solution)	7.5 - 8.5
Loss on drying	max. 0.5%
Sulphated Ash	max. 0.05%
Chloride (Cl)	max. 0.2%
Heavy Metals (Pb)	max. 0.001%

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29252990 (GST 18%)
25 Gms	29252990 (GST 18%)
500 Gms	29252990 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
5 Gms	Plastic Bottle

Available Packages

25 Gms
100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com