



19630

**L-Phenylalanine-7-Amido-4-Methylcoumarin
Trifluoroacetate Salt extrapure, 98%**

Part B

CAS : 108321-84-2

Molecular Formula : C₂₁H₁₉F₃N₂O₅

Molecular Weight : 436.4

General Information

Storage 2 to 8°C (Refrigerate)

HSN Code

25 Mg

29223900 (GST 18%)

50 Mg

29223900 (GST 18%)

Available Packages

25 Mg

50 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com