



20139

## Trimethyl Chlorosilane (TMCS) extrapure, 98%

Part B

CAS : 75-77-4

Molecular Formula : C<sub>3</sub>H<sub>9</sub>ClSi

Molecular Weight : 108.64

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay	min. 98%
Density (g/ml) @ 20°C	0.858-0.860
Refractive Index (20°C)	1.388-1.389
Boiling Range	56-58°C

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1298 , IMCO Class No.:3(8) , Packing Group:II
HSN Code	
1000 ml	29319090 (GST 18%)
25 ml	29319090 (GST 18%)
250 ml	29319090 (GST 18%)
2500 ml	29319090 (GST 18%)
Type of Packing	
1000 ml	Glass Bottle
25 ml	Glass Bottle
250 ml	Glass Bottle
2500 ml	Glass Bottle

### Available Packages

250 ml

1000 ml

2500 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com