



20770

SWCNT Type 2 Carbon Nanotubes Single-Walled (Short), Length 0.5-2.0µm, OD 1.0-2.0nm, 90%

Part C

CAS : 308068-56-6

Specifications

Appearance (Colour)	Black
Appearance (Form)	Powder
Assay	min. 90%
OD	1.0 - 2.0 nm
Length	0.5 - 2.0 µm

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1325 , IMCO Class No.:4.1 , Packing Group:III
HSN Code	
100 Mg	28030090 (GST 18%)
250 Mg	28030090 (GST 18%)

Available Packages

100 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com