



20779

## Triethanolamine pure, 98%

Part A

CAS : 102-71-6

Molecular Formula : C<sub>6</sub>H<sub>15</sub>NO<sub>3</sub>

Molecular Weight : 149.19

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to slight yellow
Appearance (Form)	Liquid
Assay (NT)	min. 98%
Density (g/ml) @ 20°C	1.122 - 1.128
Refractive Index (20°C)	1.484 - 1.488
Water (KF)	max. 0.5%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Restricted for export from India
HSN Code	
2500 ml	29221990 (GST 18%)
500 ml	29221990 (GST 18%)
Type of Packing	
2500 ml	Glass Bottle
500 ml	Glass Bottle

### Available Packages

500 ml

2500 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com