



20863

Lactophenol Mountant Solution

Part A

General Information

Storage	25 to 40°C (Room Temperature)
IMDG Identification	UN No.:2927 , IMCO Class No.:6.1(8) , Packing Group:II
HSN Code	
500 ml	32129090 (GST 18%)

Available Packages

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com