



20883

## Potassium Gluconate extrapure, 99%

Part B

CAS : 299-27-4

Molecular Formula : C<sub>6</sub>H<sub>11</sub>KO<sub>7</sub>

Molecular Weight : 234.25

### Specifications

Appearance (Colour)	White to off white
Appearance (Form)	Powder
Solubility (Turbidity) 10% aq. solution	Clear
Assay (NT) (On dry basis)	min. 99%
Loss on drying	max. 4%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29181690 (GST 18%)
500 Gms	29181690 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

### Available Packages

100 Gms  
500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com