



20975

Arachidonic Acid (non-animal) extrapure, 99%

Part B

CAS : 506-32-1

Molecular Formula : C₂₀H₃₂O₂

Molecular Weight : 304.46

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay	min. 99%

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Mg	29157090 (GST 18%)
100 Mg	29157090 (GST 18%)
25 Mg	29157090 (GST 18%)

Available Packages

10 Mg

25 Mg

100 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com