



21154

Barbituric Acid extrapure AR, 99%

Part A

CAS : 67-52-7

Molecular Formula : C₄H₄N₂O₃

Molecular Weight : 128.09

Specifications

Appearance (Colour)	White to Off-white to Cream to Light beige
Appearance (Form)	Powder
Solubility (Turbidity) 2% aq. hot solution	Clear
Assay (by Neutralization)	min. 99.0%
Melting Point	248-255°C
Loss on drying at 105°C	max. 0.5%
Sulphated Ash	max. 0.5%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29335200 (GST 18%)
25 Gms	29335200 (GST 18%)

Available Packages

25 Gms

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com