



21206

L-Lactic Acid extrapure AR, ACS, ExiPlus, Multi-Compendial, 88%

Part A

CAS : 79-33-4

Molecular Formula : C₃H₆O₃

Molecular Weight : 90.08

Specifications

ExiPlus	meets compendial specs of BP, Ph. Eur
Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Solubility with water	Miscible
Colour (APHA)	max. 10
Boiling Point	122°C
Assay (T)	min. 88%
Chloride (Cl)	max. 0.001%
Sulphate (SO ₄)	max. 0.002%
Iron (Fe)	max. 0.0001%
Heavy Metals (Pb)	max. 0.0001%
Sulphated Ash	max. 0.01%
Substance Discoloured by H ₂ SO ₄	Passes test

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
2500 ml	29181110 (GST 18%)
500 ml	29181110 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

Available Packages

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com