



21392

Methyl Acetoacetate extrapure AR, 99%

Part B

CAS : 105-45-3

Molecular Formula : C₅H₈O₃

Molecular Weight : 116.12

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 10
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	1.075-1.078
Refractive Index (20°C)	1.418-1.419
Boiling Range	168-171°C
Non Volatile Matter	max. 0.01%
Acidity (CH ₃ COOH)	max. 0.01%
Iron (Fe)	max. 0.0001%
Heavy Metals (Pb)	max. 0.0001%
Halides (Cl)	max. 0.0001%
Copper (Cu)	max. 0.00005%
Water (KF)	max. 0.05%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
2500 ml	29183040 (GST 18%)
500 ml	29183040 (GST 18%)

Available Packages

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com