



21524

Chlorophenoxyacetic Acid (4-CPA) for tissue culture, 98%

Part B

CAS : 122-88-3

Molecular Formula : C₈H₇ClO₃

Molecular Weight : 186.59

Specifications

Appearance (Colour)	White to faint beige
Appearance (Form)	Powder
Assay (GC/T)	min 98%
Loss on drying	max. 1%
Melting Point	155 - 159°C

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	98020000 (GST 18%)
25 Gms	98020000 (GST 18%)
500 Gms	98020000 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

Available Packages

25 Gms

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com