



21800

Dimethyl Dichlorosilane (DMDCS) pure, 98%

Part B

CAS : 75-78-5

Molecular Formula : C₂H₆Cl₂Si

Molecular Weight : 129.06

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay (GC)	min. 98%
Density (g/ml) @ 20°C	1.071-1.073
Refractive Index (20°C)	1.403-1.405
Boiling Range	68-70°C

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1162 , IMCO Class No.:3(8) , Packing Group:II
HSN Code	
1000 ml	29420090 (GST 18%)
250 ml	29420090 (GST 18%)
500 ml	29420090 (GST 18%)

Available Packages

250 ml

500 ml

1000 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com