



21840

Albumin Human Serum (Lipid Enhanced), 95%

Part B

CAS : 70024-90-7

Specifications

Appearance (Colour)	Pale yellow
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 5% aq. solution	Clear
Solubility (Colour) 5% aq. solution	Pale yellow
Assay (ex N)	min. 95%
pH (1% aq. solution)	6.5 - 7.5
Nitrogen (N) content	14.5% - 16.5%
Fat	min. 0.1%
Water (KF)	max. 10%

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	38229090 (GST 12%)
100 Mg	38229090 (GST 12%)
5 Gms	38229090 (GST 12%)

Available Packages

100 Mg

1 Gms

5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com