



22173

Resorcinol ACS, ExiPlus, Multi-Compendial, 99%

Part B

CAS : 108-46-3

Molecular Formula : C₆H₆O₂

Molecular Weight : 110.11

Specifications

ExiPlus	meets compendial specs of BP
Appearance (Colour)	White to slight pink
Appearance (Form)	Crystals or flakes
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Almost colourless
Assay (GC)	min. 99%
Melting Point	110 - 111°C
Sulphated Ash	max. 0.1%
Diresorcinol	max. 0.005%
Phenol	max. 0.005%
Acidity	max. 0.1ml N%
Alkalinity	max. 0.1 ml N%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2876 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
100 Gms	29072100 (GST 18%)
500 Gms	29072100 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

Available Packages

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com