



22308

Sodium Butyrate Reagent Grade, 98.5%

Part B

CAS : 156-54-7

Molecular Formula : C₄H₇O₂Na

Molecular Weight : 110.09

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline hygroscopic powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless to slight hazy
Assay (NT)	min. 98.5%
pH (2% aq. solution)	7.0 - 9.0
Water (KF)	max. 1.0%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29156010 (GST 18%)
25 Gms	29156010 (GST 18%)
500 Gms	29156010 (GST 18%)

Available Packages

25 Gms
100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com