



22350

2-Deoxyguanosine Monohydrate extrapure, 99%

Part B

CAS : 312693-72-4

Molecular Formula : C₁₀H₁₃N₅O₄.H₂O

Molecular Weight : 285.24

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Assay	max. 98%
Water (KF)	max. 7%
Spectral Ratio (A250nm/A260nm)	1.12 - 1.18
Spectral Ratio (A280nm/A260nm)	0.65 - 0.69

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29349990 (GST 18%)
100 Mg	29349990 (GST 18%)
5 Gms	29349990 (GST 18%)
500 Mg	29349990 (GST 18%)
Type of Packing	
1 Gms	Glass Bottle
100 Mg	Glass Bottle
500 Mg	Glass Bottle

Available Packages

100 Mg

500 Mg

1 Gms

5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com