



22654

2-Chloroethylphosphonic Acid (ETHREL, Ethephon) 40% soln.

Part B

CAS : 16672-87-0

Molecular Formula : C₂H₆ClO₃P

Molecular Weight : 144.50

Specifications

| | |
|------------------------|---------------------------|
| Appearance (Form) | Liquid |
| Appearance (Clarity) | Clear |
| Appearance (Colour) | Colourless to pale yellow |
| Assay (NaOH titration) | ~ 40% |

Other Information

Application

An effective germinating agent and plant growth regulator used to induce senescence and fruiting and also used as phosphorylating agent and butyrylcholinesterase inhibitor

General Information

| | |
|---------------------|--|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 60 Months |
| IMDG Identification | UN No.:2928 , IMCO Class No.:6.1(8) , Packing Group:II |
| HSN Code | |
| 100 ml | 98020000 (GST 18%) |
| 250 ml | 98020000 (GST 18%) |

Available Packages

100 ml

250 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com