



22806

Silicon Dioxide (Silica) Nanodispersion Type A (20nm), 40% in water

Part C

CAS : 7631-86-9

Molecular Formula : SiO₂

Molecular Weight : 60.08

Specifications

Appearance (Clarity)	Hazy
Appearance (Form)	Colloidal dispersion
APS	20 nm
Viscosity (20°C)	25 cPs
pH	~9.0
Concentration	40% in water

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	12 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	38229090 (GST 12%)
25 Gms	38229090 (GST 12%)
Type of Packing	
100 Gms	Plastic Bottle
25 Gms	Plastic Bottle

Available Packages

25 Gms

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com