



23576

Oil Red O for molecular biology, 75%

Part B

CAS : 1320-06-5

Molecular Formula : C₂₆H₂₄N₄O

Molecular Weight : 408.49

Specifications

Appearance (Colour)	Dark brown
Appearance (Form)	Powder
Absorbance (A) of 1% solution in Chloroform in a 1 cm cell	
@530nm	670 - 700
Loss on drying	max. 8%

Other Information

Application

Suitable for staining lipoproteins separated by cellulose acetate electrophoresis. Can be used with Agarose

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	32049000 (GST 18%)
25 Gms	32049000 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

Available Packages

25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com