



23595

Cetyltrimethyl Ammonium Chloride (CTAC)  
extrapure AR, 99%

Part B

CAS : 112-02-7

Molecular Formula : C19H42NCl

Molecular Weight : 320.00

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 5% aq. solution	Clear
Solubility (Colour) 5% aq. solution	Colourless
pH (5% aq. solution)	4.0 -7.5
Assay (NT)	min. 99%
Melting Point	232 - 237°C
Loss on drying	max.1%
Sulphate (SO4)	max. 0.01%
Iron (Fe)	max. 0.001%
Heavy Metals (Pb)	max. 0.001%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2923 , IMCO Class No.:8(6.1) , Packing Group:III
HSN Code	
100 Gms	29239000 (GST 18%)
25 Gms	29239000 (GST 18%)
500 Gms	29239000 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
25 Gms	Plastic Bottle
500 Gms	Plastic Bottle

Available Packages

25 Gms

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: [info@srlchem.com](mailto:info@srlchem.com), website [www.srlchem.com](http://www.srlchem.com)