



23622

Uracil extrapure, 99%

Part B

CAS : 66-22-8

Molecular Formula : C₄H₄N₂O₂

Molecular Weight : 112.09

Specifications

Appearance (Colour)	White to off white
Appearance (Form)	Crystalline powder
Assay (UV)	min. 99%
Loss on drying	max. 0.5%
Sulphated Ash	max. 0.1%
Heavy Metals (Pb)	max. 0.001%
Spectral Ratio (A250nm/A260nm)	0.81 - 0.85
Spectral Ratio (A280nm/A260nm)	0.18 - 0.22
Spectral Ratio (A290nm/A260nm)	max. 0.02

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kg	29333990 (GST 18%)
100 Gms	29333990 (GST 18%)
25 Gms	29333990 (GST 18%)
5 Kg	29333990 (GST 18%)
Type of Packing	
10 Gms	Plastic Bottle
100 Gms	Plastic Bottle

Available Packages

25 Gms
100 Gms
1 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com