



23700

Aluminium Chloride Hexahydrate extrapure AR,
99%

Part A

CAS : 7784-13-6

Molecular Formula : $\text{AlCl}_3 \cdot 6\text{H}_2\text{O}$

Molecular Weight : 241.43

Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline compound
Solubility (Turbidity) 10% aq. solution	Clear
Assay (Complexometric)	min. 99%
Sulphate (SO_4)	max. 0.01%
Iron (Fe)	max. 0.001%
Heavy Metals (Pb)	max. 0.001%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
500 Gms	28273200 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com