



23727

## 2-Chloropropionic Acid pure, 97%

Part B

CAS : 598-78-7

Molecular Formula : C<sub>3</sub>H<sub>5</sub>ClO<sub>2</sub>

Molecular Weight : 108.52

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay (GC)	min. 97%
Density (g/ml) @ 20°C	1.260 - 1.270
Refractive Index (20°C)	1.434 - 1.436

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2511 , IMCO Class No.:8 , Packing Group:III
HSN Code	
100 ml	29155000 (GST 18%)
500 ml	29155000 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

### Available Packages

100 ml

500 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com