



23845

## Casein acc. to Hammarsten, 95% Protein

Part B

CAS : 9000-71-9

### Specifications

Appearance (Colour)	White
Appearance (Form)	Fine powder
Assay (ex N)	min. 95% protein (dry basis)
Loss on drying	max. 6%
Sulphated Ash	max. 1%
Fat	max. 0.2%
Acidity (Lactic acid)	max. 0.1%
Alkali solubility	Passes test
Suitability as enzyme Substrate	Passes test

### Other Information

Application

Casein acc. to Hammarsten is a very high quality product for diagnostic kit production and applications. Very high solution clarity product.

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kg	35011000 (GST 18%)
100 Gms	35011000 (GST 18%)
500 Gms	35011000 (GST 18%)
Type of Packing	
1 Kg	Plastic Jar
100 Gms	Plastic Bottle

### Available Packages

100 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com