



23875

BOC-L-Amino Acids Kit

Part B

Specifications

Assay	min. 98%
Homogeneity by TLC	Passes test

Other Information

Notes

All 20 BOC-L-Amino acids, 5g each, in the kit, have assay min.99% and TLC test shows single spot.

Amino acids contained are,
BOC-L-Alanine, BOC-L-Asparagine,
I BOC-L-Arginine HCl, BOC-L-Aspartic Acid,
BOC-L-Cystine, BOC-L-Glutamic Acid,
BOC-L-Glutamine, BOC Glycine,
BOC-L-Histidine, BOC-L-Hydroxyproline,
BOC-L-Isoleucine, BOC-L-Lysine ,
BOC-L-Methionine, BOC-L-Phenylalanine,
BOC-L-Proline, BOC-L-Serine,
BOC-L-Threonine, BOC-L-Tryptophan,
BOC-L-Tyrosine, BOC-L-Valine

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kit	38229090 (GST 12%)

Available Packages

1 Kit

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com