



23882

## Tungsten Disulfide (WS2) Aerosol Spray (13Oz/369g)

Part C

CAS : 12138-09-9

### Other Information

Notes

Tungsten Disulfide (WS2) Aerosol Spray offers convenient way to apply thin film of Tungsten Disulfide (WS2) coating on any surface as well as in hard to reach places.

Tungsten Disulfide (WS2) Spray provides extreme lubricity, dry lubricant film on a substrate.

- Bonds well with substrates and form a coating film
- WS2 is chemically inert
- High Temperature lubrication
- Very high load bearing properties
- Contains no fluorocarbons or lead
- Not harmful to Ozone

The product is supplied in Aerosol Spray Cans each of 13Oz/369g.

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 piece	28309020 (GST 18%)
Type of Packing	
1 piece	Plastic Bottle

### Available Packages

1 piece

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com