



24008

Bismuth Metal Ingots ultrapure, 99.99%

Part A

CAS : 7440-69-9

Molecular Formula : Bi

Molecular Weight : 208.98

Specifications

Appearance (Form)	Shiny grey
Appearance (Colour)	Lumps
Assay	min. 99.99%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	81061020 (GST 18%)
50 Gms	81061020 (GST 18%)
Type of Packing	
50 Gms	Plastic Bottle

Available Packages

10 Gms
50 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com