



24073

Sulphamic Acid pure, 99%

Part A

CAS : 5329-14-6

Molecular Formula : H₃NO₃S

Molecular Weight : 97.09

Specifications

| | |
|------------------------|----------------------|
| Appearance (Colour) | White |
| Appearance (Form) | Crystalline compound |
| Assay (NaOH titration) | min. 99% |

General Information

| | |
|---------------------|---|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 60 Months |
| IMDG Identification | UN No.:2967 , IMCO Class No.:8 , Packing Group:III |
| HSN Code | |
| 5 Kg | 28111990 (GST 18%) |
| 500 Gms | 28111990 (GST 18%) |
| Type of Packing | |
| 5 Kg | Plastic Jar |
| 500 Gms | Plastic Bottle |

Available Packages

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com