



24167

## Boron Trioxide Anhydrous Granular (Boric Anhydride) extrapure, 97%, -60 mesh

Part A

CAS : 1303-86-2

Molecular Formula : B<sub>2</sub>O<sub>3</sub>

Molecular Weight : 69.62

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Assay	min. 97%
Iron (Fe)	max. 0.001%
Residue on ignition	max. 3%
Mesh size	min.50%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
100 Gms	28100010 (GST 18%)
250 Gms	28100010 (GST 18%)
500 Gms	28100010 (GST 18%)

### Available Packages

100 Gms  
250 Gms  
500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com