



24263

Acrylamide / Bis-acrylamide Premix Powder,
Ratio 19:1

Part A

Specifications	
Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
pH (5% aq. solution)	5.0 - 8.0
Specific conductance (25°C)	max. 10? mho/cm (5% aq. solution)
Suitability for Electrophoresis	Passes test

General Information	
Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	UN No.:2074 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
200 Gms	38229090 (GST 12%)
30 Gms	38229090 (GST 12%)
Type of Packing	
200 Gms	Glass Bottle
30 Gms	Glass Bottle

Available Packages	
30 Gms	
200 Gms	

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com