



24532

## Dichloromethane (DCM) extrapure AR, 99.5%

Part A

CAS : 75-09-2

Molecular Formula : CH<sub>2</sub>Cl<sub>2</sub>

Molecular Weight : 84.93

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 10
Assay (GC)	min. 99.5%
Density (g/ml) @ 20°C	1.323-1.325
Refractive Index (20°C)	1.424-1.425
Boiling Range	39-40°C
Non Volatile Matter	max. 0.001%
Acidity (HCl)	max. 0.001%
Free Chlorine (Cl)	max. 0.0001%
Iron (Fe)	max. 0.00005%
Heavy Metals (Pb)	max. 0.0001%
Formaldehyde (HCHO)	max. 0.0001%
Copper (Cu)	max. 0.0000 5%
Water (KF)	max. 0.05%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1593 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
25 Ltr	29031200 (GST 18%)
2500 ml	29031200 (GST 18%)
500 ml	29031200 (GST 18%)
Type of Packing	
2500 ml	Glass Bottle
500 ml	Glass Bottle

## Available Packages

500 ml

2500 ml

25 Ltr

---

### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: [info@srchem.com](mailto:info@srchem.com), website [www.srchem.com](http://www.srchem.com)

---