



24607

## Silver Iodide extrapure AR, 99%

Part A

CAS : 7783-96-2

Molecular Formula : AgI

Molecular Weight : 234.77

### Specifications

|                     |                         |
|---------------------|-------------------------|
| Appearance (Colour) | Yellowish green to grey |
| Appearance (Form)   | Powder                  |
| Assay               | min. 99%                |

### General Information

|                     |   |
|---------------------|---|
| Storage             | 25 to 40°C (Room Temperature)                         |
| Shelf Life          | 60 Months   |
| IMDG Identification | UN No.:3077 , IMCO Class No.:9 ,<br>Packing Group:III |
| HSN Code            |   |
| 10 Gms              | 28432900 (GST 18%)                                    |
| 100 Gms             | 28432900 (GST 18%)                                    |
| 25 Gms              | 28432900 (GST 18%)                                    |
| Type of Packing     |   |
| 25 Gms              | Glass Bottle  |

### Available Packages

10 Gms  
25 Gms  
100 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com