



24626

FMOC-4-Nitro-L-Phenylalanine extrapure (FMOC-L-Phe(4NO₂)-OH), 98%

Part B

CAS : 95753-55-2

Molecular Formula : C₂₄H₂₀N₂O₆

Molecular Weight : 432.43

Specifications

Appearance (Colour)	White to yellow to beige
Appearance (Form)	Flakes or powder
Assay (NT)	min. 98%

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	29379020 (GST 18%)
25 Gms	29379020 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

Available Packages

10 Gms
25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com