



24651

## Triphenylphosphine pure, 98%

Part B

CAS : 603-35-0

Molecular Formula : C<sub>18</sub>H<sub>15</sub>P

Molecular Weight : 262.30

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder/flakes
Solubility (Turbidity) 5% solution in Chloroform	Clear to slightly hazy
Solubility (Colour) 5% solution in Chloroform	Colourless
Assay (NT)	min. 98%
Melting Point	78 - 82°C

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29319090 (GST 18%)
1000 Gms	29319090 (GST 18%)
500 Gms	29319090 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
1000 Gms	Plastic Bottle

### Available Packages

100 Gms  
500 Gms  
1 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com