



24737

Tin Metal Granular extrapure AR, 99.7%

Part A

CAS : 7440-31-5

Molecular Formula : Sn

Molecular Weight : 118.69

Specifications

Appearance (Colour)	Silver grey
Appearance (Form)	Granules
Assay	min. 99.7%
Antimony (Sb)	max. 0.02%
Copper (Cu)	max. 0.005%
Iron (Fe)	max. 0.01%
Lead (Pb)	max. 0.005%
Arsenic (As)	max. 0.0001%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	80011090 (GST 18%)
250 Gms	80011090 (GST 18%)
500 Gms	80011090 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

Available Packages

100 Gms
250 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com