



25529

## Boron Trioxide Anhydrous (Boric Anhydride) extrapure AR, 98.5%

Part A

CAS : 1303-86-2

Molecular Formula : B<sub>2</sub>O<sub>3</sub>

Molecular Weight : 69.62

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Assay	min. 98.5%
Residue on ignition	max. 3%
Iron (Fe)	max. 0.001%
Sulphate (SO <sub>4</sub> )	max. 0.03%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kg	28100010 (GST 18%)
100 Gms	28100010 (GST 18%)
250 Gms	28100010 (GST 18%)

### Available Packages

100 Gms

---

250 Gms

---

1 Kg

---

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com