



25793

## TPTZ extrapure AR (2,4,6-Tri-(2-pyridyl)-5-triazine), 99%

Part B

CAS : 3682-35-7

Molecular Formula : C<sub>18</sub>H<sub>12</sub>N<sub>6</sub>

Molecular Weight : 312.34

### Specifications

Appearance (Colour)	Off - white to pale yellow
Appearance (Form)	Powder
Assay	min. 99%
Melting Point	248 - 250°C
Suitability as reagent for Iron (Fe)	Passes test

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29336990 (GST 18%)
5 Gms	29336990 (GST 18%)
Type of Packing	
1 Gms	Glass Bottle

### Available Packages

1 Gms

5 Gms

25 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com