



25818

Sodium Diphenyl Reagent AR

Part B

CAS : 5137-46-2

Molecular Formula : C₁₂H₁₀Na

Molecular Weight : 177.20

Specifications

Appearance (Colour)	Bluish black
Appearance (Form)	Viscous liquid
Activity	min. 6 meq/15ml vial
Chloride (Cl)	max. 0.0002%

Other Information

Application	a versatile reagent for rapid determination of organically bound halogens
-------------	---

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	12 Months
IMDG Identification	UN No.:2920 , IMCO Class No.:8(3) , Packing Group:II
HSN Code	
10 x15 ml	38229090 (GST 12%)
15 ml	38229090 (GST 12%)
Type of Packing	
10 x15 ml	Glass Bottle

Available Packages

15 ml
10 x15 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com