



26117

Vancomycin Supplement (3mg/vl)

Part D

Specifications

Appearance (Colour)	White, fine powder
Cultural Response	Inoculate and incubate at $35 \pm 2^\circ\text{C}$ for 24-48 hours. (Aseptically, add the contents of 1 vial of Vancomycin Supplement (26117) as per direction)
Organism	Enterococcus faecalis NCTC 12201
Inoculum (cfu) 10-100	Growth : Good
Organism	Escherichia coli ATCC 25922
Inoculum	$\geq 10^3$ Growth : Inhibited

Other Information

Applications

Used for the selective isolation of Enterococci.

Ingredients

Vancomycin	3.00mg/vial
Vancomycin	3.00mg/vial

Directions

Aseptically add the content of 1 vial to 500ml of sterile cooled at $45-50^\circ\text{C}$ VRE Agar Base or VRE Broth Base.

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 vl	38229090 (GST 12%)
Type of Packing	
5 vl	Glass Bottle

Available Packages

5 vl

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com