



26440

Petroleum Ether 40-60 extrapure AR

Part A

CAS : 8032-32-4

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Density (g/ml) @ 20°C	0.640-0.650
Boiling Range	40-60°C
Non Volatile Matter	max. 0.001%
Acidity (CH3COOH)	max. 0.001%
Iron (Fe)	max. 0.0001%
Heavy Metals (Pb)	max. 0.0001%
Bromine Number	max. 1.0
Benzene (C6H6)	max. 0.025%
Sulphur Compound (S)	max. 0.005%
Water (KF)	max. 0.01%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1268 , IMCO Class No.:3 , Packing Group:II
HSN Code	
2500 ml	27101990 (GST 18%)
500 ml	27101990 (GST 18%)
Type of Packing	
2500 ml	Glass Bottle
500 ml	Glass Bottle

Available Packages

500 ml
2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com