



26477

**N,O-Bis-(Trimethylsilyl) Acetamide (BSA)**  
**extrapure, 95%**

Part B

CAS : 10416-59-8

Molecular Formula : C<sub>8</sub>H<sub>21</sub>NOSi<sub>2</sub>

Molecular Weight : 203.43

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay (GC)	min. 95.0%

Other Information

Application

BSA is an excellent silylation reagent. For some sterically hindered primary or secondary nitro compounds, BSA and BSTFA achieve better results than normal silylation conditions.

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	UN No.:2920 , IMCO Class No.:8(3) , Packing Group:II
HSN Code	
100 ml	29319090 (GST 18%)
25 ml	29319090 (GST 18%)
250 ml	29319090 (GST 18%)

Available Packages

25 ml
100 ml
250 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com