



26500

Potassium Permanganate extrapure AR, ACS, ExiPlus, Multi-Compendial, 99.5%

Part A

CAS : 7722-64-7

Molecular Formula : $KMnO_4$

Molecular Weight : 158.03

Specifications

ExiPlus	meets compendial specs of USP, BP
Appearance (Colour)	Dark green or purple or black
Appearance (Form)	Crystals/ Granules/ Powder
Solubility (Turbidity) 5% aq. solution	Clear
Assay (Iodometric)	min. 99.5%
Insoluble matter	max. 0.2%
Chloride (Cl)	max. 0.005%
Sulphate (SO ₄)	max. 0.005%
Iron (Fe)	max. 0.002%
Heavy Metals (Pb)	max. 0.002%
Nitrogen Compounds (N)	max. 0.015%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Restricted for export from India
HSN Code	
5 Kg	28416100 (GST 18%)
500 Gms	28416100 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

Available Packages

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com