



26579

D-Tryptophan (base) extrapure, 99%

Part A

CAS : 153-94-6

Molecular Formula : C₁₁H₁₂N₂O₂

Molecular Weight : 204.23

Specifications

Appearance (Colour)	White to almost white
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% solution 2N HCl	Clear
Assay	min 99%
Specific rotation [α] _D	+31° to +32° (C=1,water)
Loss on drying	max. 0.3%
Sulphated Ash	max. 0.1%
Heavy Metals (Pb)	max. 0.001%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	29224990 (GST 18%)
25 Gms	29224990 (GST 18%)
Type of Packing	
10 Gms	Glass Bottle
25 Gms	Glass Bottle

Available Packages

10 Gms
25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com