



26645

Potassium Iodide 3M Solution

Part A

CAS : 7681-11-0

Specifications	
Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Solution
Molarity	3.0000± 0.1 M

General Information	
Storage	25 to 40°C (Room Temperature)
Shelf Life	36 Months
HSN Code	
1000 ml	38229090 (GST 12%)
500 ml	38229090 (GST 12%)

Available Packages	
500 ml	
1000 ml	

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com