



26862

O-Benzyl-L-Serine extrapure, 99%

Part B

CAS : 472696-9

Molecular Formula : C10H13NO3

Molecular Weight : 195.22

Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline compound
Assay	min. 99%
Specific rotation [?]	+20 to +22° (C=2 , 80% Acetic acid)
Melting Point (Decomposition)	~ 227°C

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29379020 (GST 18%)
10 Gms	29379020 (GST 18%)
25 Gms	29379020 (GST 18%)
5 Gms	29379020 (GST 18%)

Available Packages

1 Gms
5 Gms
10 Gms
25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com