



26892

Phenylhydrazine Hydrochloride extrapure AR, ExiPlus, Multi-Compendial, 99%

Part A

CAS : 59-88-1

Molecular Formula : C₆H₈N₂.HCl

Molecular Weight : 144.60

Specifications

ExiPlus	meets compendial specs of USP, BP, Ph.Eur
Appearance (Colour)	White to off - white
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 2% aq. solution	Clear
Assay (T)	min. 99%
Melting Point	243 - 246°C
Sulphated Ash	max. 0.05%
Sulphate (SO ₄)	max. 0.01%
Iron (Fe)	max. 0.001%
Heavy Metals (Pb)	max. 0.001%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2811 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
100 Gms	29280090 (GST 18%)
250 Gms	29280090 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
250 Gms	Glass Bottle

Available Packages

100 Gms

250 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com