



27166

3-Nitrophenylboronic Acid extrapure, 97%

Part B

CAS : 13331-27-6

Molecular Formula : C₆H₆BN₂O₄

Molecular Weight : 166.9

Specifications

Appearance (Colour)	White to pale yellow
Appearance (Form)	Crystalline powder
Assay (NaOH titration)	min. 97%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29319090 (GST 18%)
25 Gms	29319090 (GST 18%)
5 Gms	29319090 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

Available Packages

1 Gms

5 Gms

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com