



27297

Bathocuproine (BCP) ACS, 99%

Part B

CAS : 4733-39-5

Molecular Formula : C₂₆H₂₀N₂

Molecular Weight : 360.46

Specifications

Appearance (Colour)	Off white to buff
Appearance (Form)	Powder
Assay (NT)	min. 99%
Melting Point	279 - 283°C
Clarity of solution	Passes test
Sensitivity to copper (Cu)	Passes test

Other Information

Application	Reagent for determination of Copper (Cu).
-------------	---

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
100 Mg	29335990 (GST 18%)
250 Mg	29335990 (GST 18%)

Available Packages

500 Mg
1 Gms
5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com